

Product datasheet for **RC209399L4V**

GTF2H3 (NM_001516) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GTF2H3 (NM_001516) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GTF2H3
Synonyms:	BTF2; P34; TFB4; TFIIH
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001516
ORF Size:	924 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209399).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001516.3
RefSeq Size:	3434 bp
RefSeq ORF:	927 bp
Locus ID:	2967
UniProt ID:	Q13889
Cytogenetics:	12q24.31
Domains:	Tfb4
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

Protein Pathways: Basal transcription factors, Nucleotide excision repair

MW: 34.4 kDa

Gene Summary: This gene encodes a member of the TFB4 family. The encoded protein is a subunit of the core-TFIID basal transcription factor and localizes to the nucleus. The encoded protein is involved in RNA transcription by RNA polymerase II and nucleotide excision repair and associates with the Cdk-activating kinase complex. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 14. [provided by RefSeq, Dec 2012]