

Product datasheet for **RC209426L3V**

SSBP3 (NM_145716) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SSBP3 (NM_145716) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SSBP3
Synonyms:	CSDP; SSDP; SSDP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_145716
ORF Size:	1164 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209426).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145716.2
RefSeq Size:	3259 bp
RefSeq ORF:	1167 bp
Locus ID:	23648
UniProt ID:	Q9BWW4
Cytogenetics:	1p32.3
Domains:	LisH, SSDP
Protein Families:	Transcription Factors



[View online »](#)

MW: 40.4 kDa

Gene Summary: May be involved in transcription regulation of the alpha 2(I) collagen gene where it binds to the single-stranded polypyrimidine sequences in the promoter region.[UniProtKB/Swiss-Prot Function]