

Product datasheet for **RC209450A1V**

Human PAF Receptor (PTAFR) (NM_000952) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human PAF Receptor (PTAFR) (NM_000952) AAV Particle
Tag:	Myc-DDK
Symbol:	PAF Receptor
Synonyms:	PAFR
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC209450 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGCCACATGACTCCTCCACATGGACTCTGAGTTCGGATACACTCTCTCCCGATTGTTTACAGCA
TCATCTTTGTGCTCGGGTCATTGCTAATGGCTACGTGCTGTGGGTCTTTGCCCGCTGTACCCTTGCAA
GAAATTCATGAGATAAAGATCTTCATGGTGAACCTCACCATGGCGGACATGCTCTTCTGATCACCCCTG
CCACTTTGGATTGTCTACTACAAAACAGGGCAACTGGATACTCCCAAATTCCTGTGCAACGTGGCTG
GCTGCCTTTTCTCATCAACACCTACTGCTCTGTGGCCTCCTGGCGTCACTACTATAACCGTTCCA
GGCAGTAACTCGGCCATCAAGACTGCTCAGGCCAACCCGCAAGCGTGGCATCTCTTTGTCTTGTC
ATCTGGGTGGCCATTGTGGGAGCTGCATCCTACTTCTCATCCTGGACTCCACCAACACAGTGCCCGACA
GTGCTGGCTCAGGCAACGTCACTCGCTGCTTTGAGCATTACGAGAAGGGCAGCGTGCCAGTCTCATCAT
CCACATCTTCATCGTGTTCAGCTTCTTCTGGTCTTCTCATCATCCTCTTCTGCAACCTGGTCAATC
CGTACCTTGCTCATGCAGCCGGTGCAGCAGCAGCGCAACGCTGAAGTCAAGCGCCGGGCGCTGTGGATG
TGTGCACGGTCTTGGCGGTGTTTCATCATCTGCTTCGTGCCACCACGTGGTGCAGTGCCTGGACCCT
TGCTGAGCTGGGCTTCCAGGACAGCAAATTCACCAGGCCATTAATGATGCACATCAGGTACCCTCTGC
CTCCTTAGCACCAACTGTGTCTTAGACCCTGTTATCTACTGTTTCTCACCAGAAGTTCGCAAGCACC
TCACCGAAAAGTTCTACAGCATGCGCAGTAGCCGAAAATGCTCCCGGGCCACCACGGATACGGTCACTGA
AGTGGTTGTGCCATTCAACCAGATCCCTGGCAATTCCTCAAAAAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTAA



[View online »](#)

Protein Sequence:	<p>>RC209450 protein sequence Red=Cloning site Green=Tags(s)</p> <p>MEPHDSSSHMDSEFRYTLFPVIVYSIIFVLGVIANGYVLWVWFARLYPCKKFNEIKIFMVNLTADMLFLITL PLWIVYYQNQGNWILPKFLCNVAGCLFFINTYCSVAFLGVITYNRFQAVTRPIKTAQANTRKRGISLSLV IWWAIVGAASYFLILDSTNTVPDSAGSGNVTRCFEYHEKGSVPVLIHIFIVFSFLLVFLIILFCNLVII RTLLMQPVQQRNAEVKRRALWMVCTVLAVFIICFVPHHVVQLPWTLAELGFQDSKFHQAINDAHQVTLCLLSTNCVLD PVIYCFLLTKKFRKHLTEKFYSMRSSRKCSRATTDTVTEVVVFPNQIPGNSLKN</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Species:	Human
Serotype:	AAV-2
ACCN:	NM_000952
ORF Size:	1026 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_000952.3
RefSeq Size:	4247 bp
RefSeq ORF:	1029 bp
Locus ID:	5724
UniProt ID:	P25105
Cytogenetics:	1p35.3
MW:	39.2 kDa