

Product datasheet for **RC209455L3V**

Cytochrome b reductase 1 (CYBRD1) (NM_024843) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cytochrome b reductase 1 (CYBRD1) (NM_024843) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Cytochrome b reductase 1
Synonyms:	CYB561A2; DCYTB; FRRS3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024843
ORF Size:	858 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209455).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024843.2
RefSeq Size:	4390 bp
RefSeq ORF:	861 bp
Locus ID:	79901
UniProt ID:	Q53TN4
Cytogenetics:	2q31.1
Domains:	Cytochrome_B561, B561



[View online »](#)

Protein Families:

Transmembrane

MW:

31.7 kDa

Gene Summary:

This gene is a member of the cytochrome b(561) family that encodes an iron-regulated protein. It highly expressed in the duodenal brush border membrane. It has ferric reductase activity and is believed to play a physiological role in dietary iron absorption. [provided by RefSeq, Jul 2008]