

Product datasheet for **RC209490A1V**

Human C1orf198 (NM_032800) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human C1orf198 (NM_032800) AAV Particle
Tag: Myc-DDK
Symbol: C1orf198
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC209490 representing NM_032800
Red=Cloning site Blue=ORF Green=Tags(s)

CTATAGGGCGGCCGGAATTTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCCGCGCGC

ATGGCGTCCATGGCGGCGGCGATCGCGGCTTCGCGCTCGGCGGTATGAGCGGGAACCGGCCTCTGGACG
ACCGGGAGCGAAAGCGCTTCACTTACTTCTCGTCGCTGAGCCCCATGGCCAGGAAGATCATGCAGGACAA
GGAGAAGATCCGCGAGAAGTACGGGCCGAGTGGGCGCGGCTGCCGCCGCGCAGCAGGACGAGATCATC
GACCGGTGCCTGGTGGGCGCGCGCCCCGGCGCCCCGAGACCCGGGGACTCGGAGGAGCTCACGCGCT
TCCCCGGCTTGC GCGGGCCACGGGCCAGAAGTGGTGCCTTCGGGACGAGGATCTAACTTGGCAAGA
TGAGCACTCTGCCCTTTCTCCTGGGAAACAAAGATCAGATGGAGTTTCACTATCTCCGCCCTATCCATC
CAGGAGCCGAGCAACGGCACCGCCGCGCAGCAGCCAGACCACTGTCCAAAGTTCCAGGGCTCCCAAGG
CCCTCAAGTCCTCCCAAGGCAGCAGGTCTCCAGCCTGGACGCCCTGGGCCCCACCAGGAAGGAGGAGGA
AGCGTCATTCTGGAAGATCAATGCTGAGCGGTCCCGAGGGGAGGGGCTGAGGCCGAGTTCCAGTCGCTG
ACCCCTAGCCAGATCAAGTCCATGGAGAAGGGGAAAAGGTCTTGCTCCCTGCTACCGGCAGGAACCTG
CCCCGAAGGACAGGGAGGCCAAGGTGAAAGGCCAGCACCTCCGTGAGGAGCAGCGTCCTCTTCCAA
CGTGAGCACCGAACGTGAGAGACCCAGCCTGTCCAGGCCTTCAGCAGTGCACGAGGCTGCCCCC
TCCAGCTCGAGGGGAAGCTGCCATCTCCTGATGTGAGCAGGACGATGGGAAGACACCTGTCTCGG
AACCCAAGTTTGACAGGTGAGTCAAGTAATGTCGCTTGAAGACGGGATTTGATTTTCTGGACAATTG
G

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC209490 representing NM_032800 Red=Cloning site Green=Tags(s) MASMAAAIAASRSVMSGNRPLDDRERKRFTYFSSLSMARKIMQDKEKIREKYGPEWARLPPAQQDEII DRCLVGPRAPAPRDPGDSEELTRFPGLRGPTGQKVVRFGDEDLTWQDEHSAPFSWETKSQMEFSISALSI QEPSNGTAASEPRPLSKASQGSQALKSSQGSRSSLDALGPTRKEEEASFVKINAERSRGEGPEAEFQSL TPSQIKSMEKGEKVLPPCYRQEPAPKDREKVERPSTLRQEQRPLPNVSTERERPQPVQAFSSALHEAAP SQLEGLKLPSPDVRQDDGEDTLFSEPKFAQVSSSNVVLKTGFDFLDNW TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_032800
ORF Size:	981 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_032800.1</u>
RefSeq Size:	3733 bp
RefSeq ORF:	984 bp
Locus ID:	84886
UniProt ID:	<u>Q9H425</u>
Cytogenetics:	1q42.2
MW:	36.2 kDa