

## Product datasheet for **RC209535L4V**

### **TOMM70A (TOMM70) (NM\_014820) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	TOMM70A (TOMM70) (NM_014820) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TOMM70A
Synonyms:	Tom70; TOMM70A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014820
ORF Size:	1824 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209535).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_014820.4</a>
RefSeq Size:	4381 bp
RefSeq ORF:	1827 bp
Locus ID:	9868
UniProt ID:	<a href="#">O94826</a>
Cytogenetics:	3q12.2
Domains:	TPR
Protein Families:	Transmembrane



[View online »](#)

**MW:** 67.5 kDa

**Gene Summary:** This gene encodes an import receptor of the outer mitochondrial membrane that is part of the translocase of the outer membrane complex. This protein is involved in the import of mitochondrial precursor proteins. [provided by RefSeq, Oct 2011]