

Product datasheet for RC209575L4V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Amyloid Precursor Protein (APP) (NM 201413) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: Amyloid Precursor Protein (APP) (NM 201413) Human Tagged ORF Clone Lentiviral Particle

Symbol: Amyloid Precursor Protein

Synonyms: AAA; ABETA; ABPP; AD1; alpha-sAPP; APPI; CTFgamma; CVAP; PN-II; PN2; preA4

Mammalian Cell

Selection:

Puromycin

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

Tag: mGFP

ACCN: NM_201413 **ORF Size:** 2253 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC209575).

Sequence:
OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeq: <u>NM 201413.1</u>

RefSeq Size:3584 bpRefSeq ORF:2256 bp

Locus ID: 351

UniProt ID: P05067
Cytogenetics: 21q21.3

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Alzheimer's disease





Amyloid Precursor Protein (APP) (NM_201413) Human Tagged ORF Clone Lentiviral Particle – RC209575L4V

MW: 83 kDa

Gene Summary:

This gene encodes a cell surface receptor and transmembrane precursor protein that is cleaved by secretases to form a number of peptides. Some of these peptides are secreted and can bind to the acetyltransferase complex APBB1/TIP60 to promote transcriptional activation, while others form the protein basis of the amyloid plaques found in the brains of patients with Alzheimer disease. In addition, two of the peptides are antimicrobial peptides, having been shown to have bacteriocidal and antifungal activities. Mutations in this gene have been implicated in autosomal dominant Alzheimer disease and cerebroarterial amyloidosis (cerebral amyloid angiopathy). Multiple transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Aug 2014]