

Product datasheet for **RC209588L4V**

OLFML2B (NM_015441) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	OLFML2B (NM_015441) Human Tagged ORF Clone Lentiviral Particle
Symbol:	OLFML2B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_015441
ORF Size:	2250 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209588).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015441.1 , NP_056256.1
RefSeq Size:	3796 bp
RefSeq ORF:	2253 bp
Locus ID:	25903
UniProt ID:	Q68BL8
Cytogenetics:	1q23.3
Protein Families:	Druggable Genome, Transmembrane
MW:	84 kDa


[View online »](#)

Gene Summary:

This gene encodes an olfactomedin domain-containing protein. Most olfactomedin domain-containing proteins are secreted glycoproteins. [provided by RefSeq, Dec 2016]