

Product datasheet for **RC209589A1V**

Human ZNF784 (NM_203374) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ZNF784 (NM_203374) AAV Particle
Tag:	Myc-DDK
Symbol:	ZNF784
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC209589 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCCGCTGCGCGCCAGAGGCCAGAGTCCGAGCTCACCAGCTCCGGAGTCGCGATCCCAGGAGCCAC
TGGACCTGGTCCTGGTGCCTGATGACTGCCGGCTGGCACACCCCGAGTGACCTCATCGAGATCCAGGT
GGTGAAGGTGACGGACACCAGCTGGTCCCTGAGCCCCGGAGCCAGGTTCTTTCCACTGTGCCTTGTGC
CCTGCTGCCTTCCGGCTGGTTTCCGAGCTGCTGTTCCACGAACATGGCCACTTGGCTGGGGCCGAGGGAG
GCGGGCAGGGTGGGGACCCGAGCCGGTGCCACGTGTGCGGCCACAGCTGCCCGGGCCCCGCCAGCCTGCG
CGCGCACTACAGCTTGACACGGGGGAGCGGCCCTACCGCTGCGCGCTCTGCCCGCGCGCTTCAAGGCC
TTGGCGCCCCCTGCTCCGGCACCAGCACCGGCACGGGTGGAGCCGGGCACCTCTCGGAGGCCCTCCGACA
CAGCGGCCGTTGAGAACAGAGGCCCGGGTGGCCCCGGAGAGGGCGGAGGTGGTGATGGCGGCGCGCGC
GGCGGGCGCAGCGGTGGGAAGCCTTTTGCCTGCAGGTTCTGCGCAAGCCCTTCCGCCGCTCCTCAGAC
ATGCGAGACCACGAGCGCTGCACACTGGCGAGCGGCCATACCATGCGGCATCTGCGCAAGGGCTTCA
CCCAGTCCTCGGTGCTTAGCGGCCACGCCCGCATCCACACTGGCGAGCGCCGTTCCGCTGCACGCTCTG
CGACCGCACCTTCAACAACCTCTCCAATTCCGAAGCACCAGCGCACCACTTCCACGGGCCGGGGCCG
GGGCTGGGAGACTCTGGAGGCCAGCTGGGCTCGTCGGCGGCTGAGGGGTGCGGGAGCGGGTGTGGGTAG
GGGACCTGCGGAGGAGGGGCGGGGGAGACCGCAAGGTGAAGGTGGAGGCCGACCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC209589 protein sequence Red=Cloning site Green=Tags(s) MAAARPEAQRSSPTPESRSQEPLDLVLPDDCRPGTPPSDLIEIQVVKVTDTTLVPEPPEPGSFHCALC PAAFRLVSELLFHEHGLAGAEGGGQGGDPSRCHVCGHSCPGPASLRAHYSLHTGERPYRCALCPRAFKA LAPLLRHQHRHGVPGTSRRPPDTAAVAEQRPVAPERAEEVMAAAAAGAAVGKPFACRFCAKPFRRSSD MRDHERVHTGERPYHCGICGKGFTQSSVLSGHARIHTGERPFRCTLCDRTFNSSNFRKHQRTHFHGPGP GLGDSGGQLGSSAAEGSGSGCGVGDPAEEGRGETAKVKVEADQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_203374
ORF Size:	969 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_203374.1</u>
RefSeq Size:	2006 bp
RefSeq ORF:	972 bp
Locus ID:	147808
UniProt ID:	<u>Q8NCA9</u>
Cytogenetics:	19q13.42
MW:	34.2 kDa