

## Product datasheet for **RC209601L4V**

### TRIM9 (NM\_015163) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TRIM9 (NM_015163) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TRIM9
Synonyms:	RNF91; SPRING
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_015163
ORF Size:	2130 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209601).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_015163.4</a> , <a href="#">NP_055978.3</a>
RefSeq Size:	4946 bp
RefSeq ORF:	2133 bp
Locus ID:	114088
UniProt ID:	<a href="#">Q9C026</a>
Cytogenetics:	14q22.1
Domains:	zf-B_box, BBC, SPRY, FN3
Protein Families:	Druggable Genome



[View online »](#)

**MW:** 79.2 kDa

**Gene Summary:** The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The protein localizes to cytoplasmic bodies. Its function has not been identified. Alternate splicing of this gene generates two transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]