

Product datasheet for **RC209618L3V**

ADPRH (NM_001125) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ADPRH (NM_001125) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ADPRH
Synonyms:	ARH1; hARH1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001125
ORF Size:	1071 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209618).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001125.2
RefSeq Size:	3704 bp
RefSeq ORF:	1074 bp
Locus ID:	141
UniProt ID:	P54922
Cytogenetics:	3q13.33
Domains:	ADP-ribosyl_GH
MW:	39.5 kDa



[View online »](#)

Gene Summary:

The enzyme encoded by this gene catalyzes removal of mono-ADP-ribose from arginine residues of proteins in the ADP-ribosylation cycle. Unlike the rat and mouse enzymes that require DTT for maximal activity, the human enzyme is DTT-independent. Alternatively spliced transcript variants that encode different protein isoforms have been described. [provided by RefSeq, May 2014]