

Product datasheet for **RC209630L3V**

TRIM11 (NM_145214) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TRIM11 (NM_145214) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TRIM11
Synonyms:	BIA1; RNF92
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_145214
ORF Size:	1404 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209630).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145214.2
RefSeq Size:	2715 bp
RefSeq ORF:	1407 bp
Locus ID:	81559
UniProt ID:	Q96F44
Cytogenetics:	1q42.13
Domains:	zf-B_box, RING, SPRY, PRY
Protein Families:	Druggable Genome


[View online »](#)

MW: 52.8 kDa

Gene Summary: The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. This protein localizes to the nucleus and the cytoplasm. Its function has not been identified. [provided by RefSeq, Jul 2008]