

## Product datasheet for **RC209655A1V**

### Human ENSA (NM\_207168) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human ENSA (NM\_207168) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** ENSA  
**Synonyms:** ARPP-19e  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC209655 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTCCAGAAACAAGAAGAAGAGAACCCTGCGGAGGAGACCGGCGAGGAGAAGCAGGACACGCAGGAGA  
AAGAAGGTATTCTGCCTGAGAGAGCTGAAGAGGCAAAGCTAAAGGCCAAATACCCAAGCCTAGGACAAAA  
GCCTGGAGGCTCCGACTTCTCATGAAGAGACTCCAGAAAGGGGTATGGGGCATAGTCTCTTACCCTCTT  
TCTTTGGAGCTAAAGGAGGTTCTCGAATGAAGTCTGTAGAGGTTCTACTGGATCCTTTCTGGAGGTTT  
TGTTGTTGAACAGAAGTAGAGGCGAATTTGAAATT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC209655 protein sequence  
Red=Cloning site Green=Tags(s)  
MSQKQEEENPAEETGEEKQDTQEKEGILPERAEAKLKAKYPSLGQKPGGSDFLMKRLQKGVWGI VSYPL  
SLELKEVLRMKSVEVLLDPFLEVLLLNRSRGEFEI

**TR**TRPLEQKLI SEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_207168  
**ORF Size:** 315 bp



[View online »](#)

<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u>NM_207168.1, NP_997051.1</u>
<b>RefSeq Size:</b>	771 bp
<b>RefSeq ORF:</b>	318 bp
<b>Locus ID:</b>	2029
<b>UniProt ID:</b>	<u>O43768</u>
<b>Cytogenetics:</b>	1q21.3
<b>MW:</b>	12 kDa