

Product datasheet for **RC209667**

SFT2B (SFT2D2) (NM_199344) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: SFT2B (SFT2D2) (NM_199344) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: SFT2B
Synonyms: dj747L4.C1.2; UNQ512
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC209667 representing NM_199344
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACAAGCTGAAGAAGGTGCTGAGCGGGCAGGACACGGAGGACCGGAGCGGCCTGCCAGGTTGTTG
AGGCATCTTCATTAAGCTGGAGTACCAGGATAAAAGGCTTCATTGCGTGTTCATAGGAATTCCTG
CTCACTGCTGGTACTGTTCTGCTGTGGTGCCAGGAAGGGACTACACCTCTTCGCAGTGTTCATACC
TTTGTAATATCGCATCAATTGGGAGTACCATCTCCTCATGGGACCAGTAAACAGCTGAAGCGAATGT
TTGAGCCTACTCGTTTGATTGCAACTATCATGGTGTGTTGTGTTTGCACCTACCCTGTGTTCTGCCTT
TTGGTGGCATAACAAGGGACTTGCACTTATCTTCTGCATTTTGCAGTCTTTGGCATTGACGTGGTACAGC
CTTTCCTTCATACCATTGCAAGGGATGCTGTGAAGAAGTGTTCCTGCTGTCTTGCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC209667 representing NM_199344
Red=Cloning site Green=Tags(s)

MDLKKVLSGQDTEDRSGLSEVVEASSLSWSTRIKGFIACFIAIGILCSLLGTVLLWVPRKGLHLFAVFYT
FGNIASIGSTIFLMGPVKQLKRMFEPTRLIATIMVLLCFALTLCSAFWWHNKGLALIFCQLSLALTWYS
LSFIPFARDAVKKCFVCLA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mg3781_e05.zip



[View online >](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_199344

ORF Size: 480 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_199344.3](#)


RefSeq Size: 1155 bp

RefSeq ORF: 483 bp


Locus ID: 375035

UniProt ID: [O95562](#)
 Cytogenetics: 1q24.2
 Protein Families: Transmembrane
 MW: 17.6 kDa
 Gene Summary: May be involved in fusion of retrograde transport vesicles derived from an endocytic compartment with the Golgi complex.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC209667



Western blot validation of overexpression lysate (Cat# [LY404531]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC209667 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).