

Product datasheet for **RC209667L4V**

SFT2B (SFT2D2) (NM_199344) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SFT2B (SFT2D2) (NM_199344) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SFT2B
Synonyms:	dj747L4.C1.2; UNQ512
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_199344
ORF Size:	480 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209667).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_199344.1
RefSeq Size:	1155 bp
RefSeq ORF:	483 bp
Locus ID:	375035
UniProt ID:	O95562
Cytogenetics:	1q24.2
Protein Families:	Transmembrane
MW:	17.6 kDa



[View online »](#)

Gene Summary:

May be involved in fusion of retrograde transport vesicles derived from an endocytic compartment with the Golgi complex.[UniProtKB/Swiss-Prot Function]