

Product datasheet for **RC209674L1V**

GPR83 (NM_016540) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GPR83 (NM_016540) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GPR83
Synonyms:	GIR; GPR72
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_016540
ORF Size:	1269 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209674).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016540.2
RefSeq Size:	4300 bp
RefSeq ORF:	1272 bp
Locus ID:	10888
UniProt ID:	Q9NYM4
Cytogenetics:	11q21
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Neuroactive ligand-receptor interaction



[View online »](#)

MW: 48.3 kDa

Gene Summary: Orphan receptor. Could be a neuropeptide Y receptor.[UniProtKB/Swiss-Prot Function]