

Product datasheet for **RC209696L1V**

RPS27A (NM_002954) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | RPS27A (NM_002954) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | RPS27A |
| Synonyms: | CEP80; HEL112; S27A; UBA80; UBC; UBCEP1; UBCEP80 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_002954 |
| ORF Size: | 468 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC209696). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_002954.3 |
| RefSeq Size: | 1068 bp |
| RefSeq ORF: | 471 bp |
| Locus ID: | 6233 |
| UniProt ID: | P62979 |
| Cytogenetics: | 2p16.1 |
| Domains: | UBQ, Ribosomal_S27 |
| Protein Families: | Druggable Genome |



[View online »](#)

Protein Pathways: Ribosome

MW: 18 kDa

Gene Summary: Ubiquitin, a highly conserved protein that has a major role in targeting cellular proteins for degradation by the 26S proteasome, is synthesized as a precursor protein consisting of either polyubiquitin chains or a single ubiquitin fused to an unrelated protein. This gene encodes a fusion protein consisting of ubiquitin at the N terminus and ribosomal protein S27a at the C terminus. When expressed in yeast, the protein is post-translationally processed, generating free ubiquitin monomer and ribosomal protein S27a. Ribosomal protein S27a is a component of the 40S subunit of the ribosome and belongs to the S27AE family of ribosomal proteins. It contains C4-type zinc finger domains and is located in the cytoplasm. Pseudogenes derived from this gene are present in the genome. As with ribosomal protein S27a, ribosomal protein L40 is also synthesized as a fusion protein with ubiquitin; similarly, ribosomal protein S30 is synthesized as a fusion protein with the ubiquitin-like protein fubi. Multiple alternatively spliced transcript variants that encode the same proteins have been identified.[provided by RefSeq, Sep 2008]