

Product datasheet for **RC209707L2V**

Cytokeratin 19 (KRT19) (NM_002276) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Cytokeratin 19 (KRT19) (NM_002276) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | Cytokeratin 19 |
| Synonyms: | CK19; K1CS; K19 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-mGFP (PS100071) |
| Tag: | mGFP |
| ACCN: | NM_002276 |
| ORF Size: | 1200 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC209707). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_002276.4 |
| RefSeq Size: | 1490 bp |
| RefSeq ORF: | 1203 bp |
| Locus ID: | 3880 |
| UniProt ID: | P08727 |
| Cytogenetics: | 17q21.2 |
| Domains: | filament |
| MW: | 43.9 kDa |



[View online »](#)

Gene Summary:

The protein encoded by this gene is a member of the keratin family. The keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into cytokeratins and hair keratins. The type I cytokeratins consist of acidic proteins which are arranged in pairs of heterotypic keratin chains. Unlike its related family members, this smallest known acidic cytokeratin is not paired with a basic cytokeratin in epithelial cells. It is specifically expressed in the periderm, the transiently superficial layer that envelopes the developing epidermis. The type I cytokeratins are clustered in a region of chromosome 17q12-q21. [provided by RefSeq, Jul 2008]