

Product datasheet for **RC209719A1V**

Human **CDC26** (NM_139286) AAV Particle

Product data:

Product Type:	AAV Particles
Tag:	Myc-DDK
Symbol:	CDC26
Synonyms:	ANAPC12; APC12; C9orf17
Mammalian Cell	None
Selection:	
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC209719 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s) TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC GCCCGCATCGCC ATGCTGAGACGGAAACCAACACGCCTAGAGCTAAAGCTTGATGACATTGAAGAGTTTGAGAACATTCGAA AGGACCTGGAGACCCGTAAGAAACAGAAGGAAGATGTGGAAGTTGTAGGAGGCAGTGATGGAGAAGGAGC CATTGGGCTTAGCAGTGATCCCAAGAGCCGGGAACAAATGATCAATGATCGGATTGGTTATAAACCCCAA CCCAAGCCCAATAATCGTTCATCTCAATTTGGAAGTCTTGAATTT ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT ACAAGGATGACGACGATAAGGTTTAA
Protein Sequence:	>RC209719 protein sequence Red=Cloning site Green=Tags(s) MLRRKPTRLKLDDIEEFENIRKDLETRKKQKEDVEVVGSDGEGAIGLSSDPKSREQMINDRIGYKPQ PKPNNRSSQFGSLEF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_139286
ORF Size:	255 bp



Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_139286.3</u>
RefSeq Size:	885 bp
RefSeq ORF:	258 bp
Locus ID:	246184
UniProt ID:	<u>Q8NHZ8</u>
Cytogenetics:	9q32
MW:	9.8 kDa