

## Product datasheet for **RC209748L2V**

### PRPK (TP53RK) (NM\_033550) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PRPK (TP53RK) (NM_033550) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PRPK
Synonyms:	BUD32; C20orf64; dj101A2; GAMOS4; Nori-2; Nori-2p; PRPK; TPRKB
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_033550
ORF Size:	759 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209748).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_033550.3</a>
RefSeq Size:	3384 bp
RefSeq ORF:	762 bp
Locus ID:	112858
UniProt ID:	<a href="#">Q96S44</a>
Cytogenetics:	20q13.12
Domains:	S_TKc, KOW
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

MW: 28.1 kDa

**Gene Summary:** Component of the EKC/KEOPS complex that is required for the formation of a threonylcarbamoyl group on adenosine at position 37 (t(6)A37) in tRNAs that read codons beginning with adenine (PubMed:22912744, PubMed:27903914). The complex is probably involved in the transfer of the threonylcarbamoyl moiety of threonylcarbamoyl-AMP (TC-AMP) to the N6 group of A37 (PubMed:22912744, PubMed:27903914). TP53RK has ATPase activity in the context of the EKC/KEOPS complex and likely plays a supporting role to the catalytic subunit OSGEP (By similarity). Atypical protein kinase that phosphorylates 'Ser-15' of p53/TP53 protein and may therefore participate in its activation (PubMed:11546806). [UniProtKB/Swiss-Prot Function]