

Product datasheet for RC209761L3V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

VCAM1 (NM 001078) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: VCAM1 (NM_001078) Human Tagged ORF Clone Lentiviral Particle

Symbol:

CD106: INCAM-100 Synonyms:

Mammalian Cell

Selection:

Puromycin

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

Tag: Myc-DDK NM 001078 ACCN: **ORF Size:** 2217 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC209761).

OTI Disclaimer:

Sequence:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeq: NM 001078.2, NP 001069.1

RefSeq Size: 3220 bp RefSeq ORF: 2220 bp Locus ID: 7412 **UniProt ID:** P19320 Cytogenetics: 1p21.2

Domains: ig, IGv, IGc2, IG

Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Transmembrane





Protein Pathways: Cell adhesion molecules (CAMs), Leukocyte transendothelial migration


MW: 81.3 kDa

Gene Summary: This gene is a member of the Ig superfamily and encodes a cell surface sialoglycoprotein

expressed by cytokine-activated endothelium. This type I membrane protein mediates leukocyte-endothelial cell adhesion and signal transduction, and may play a role in the development of artherosclerosis and rheumatoid arthritis. Three alternatively spliced transcripts encoding different isoforms have been described for this gene. [provided by

RefSeq, Dec 2010]

Product images:



[RC209761L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC209761L3V particle to overexpress human VCAM1-Myc-DDK fusion protein.