

Product datasheet for **RC209763L3V**

PP4R4 (PPP4R4) (NM_058237) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PP4R4 (PPP4R4) (NM_058237) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PP4R4
Synonyms:	CFAP14; KIAA1622; PP4R4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_058237
ORF Size:	2619 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209763).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_058237.1
RefSeq Size:	3868 bp
RefSeq ORF:	2622 bp
Locus ID:	57718
UniProt ID:	Q6NUP7
Cytogenetics:	14q32.12-q32.13
Domains:	HEAT
Protein Families:	Phosphatase



[View online »](#)

MW: 99.5 kDa

Gene Summary: The protein encoded by this gene is a HEAT-like repeat-containing protein. The HEAT repeat is a tandemly repeated, 37-47 amino acid long module occurring in a number of cytoplasmic proteins. Arrays of HEAT repeats form a rod-like helical structure and appear to function as protein-protein interaction surfaces. The repeat-containing region of this protein has some similarity to the constant regulatory domain of the protein phosphatase 2A PR65/A subunit. The encoded protein binds protein serine/threonine phosphatase 4c in the cytoplasm. [provided by RefSeq, Jan 2017]