

Product datasheet for **RC209775L2V**

RPS13 (NM_001017) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | RPS13 (NM_001017) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | RPS13 |
| Synonyms: | S13 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-mGFP (PS100071) |
| Tag: | mGFP |
| ACCN: | NM_001017 |
| ORF Size: | 453 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC209775). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001017.2 |
| RefSeq Size: | 529 bp |
| RefSeq ORF: | 456 bp |
| Locus ID: | 6207 |
| UniProt ID: | P62277 |
| Cytogenetics: | 11p15.1 |
| Domains: | Ribosomal_S15 |
| Protein Pathways: | Ribosome |



[View online »](#)

MW: 17 kDa

Gene Summary: Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 40S subunit. The protein belongs to the S15P family of ribosomal proteins. It is located in the cytoplasm. The protein has been shown to bind to the 5.8S rRNA in rat. The gene product of the E. coli ortholog (ribosomal protein S15) functions at early steps in ribosome assembly. This gene is co-transcribed with two U14 small nucleolar RNA genes, which are located in its third and fifth introns. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. [provided by RefSeq, Jul 2008]