

Product datasheet for **RC209783L4V**

C10orf97 (FAM188A) (NM_024948) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C10orf97 (FAM188A) (NM_024948) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C10orf97
Synonyms:	C10orf97; CARP; DERP5; FAM188A; MST126; MSTP126; my042
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024948
ORF Size:	1335 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209783).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024948.2
RefSeq Size:	2375 bp
RefSeq ORF:	1338 bp
Locus ID:	80013
UniProt ID:	Q9H8M7
Cytogenetics:	10p13
Protein Families:	Druggable Genome
MW:	49.7 kDa


[View online »](#)

Gene Summary:

The protein encoded by this gene contains a caspase-associated recruitment domain and may function in apoptosis. It has been identified as a tumor suppressor in lung and gastric cancers, and a polymorphism in the gene may be associated with gastric cancer risk. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Dec 2015]