

## Product datasheet for **RC209792L4V**

### DDX43 (NM\_018665) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	DDX43 (NM_018665) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DDX43
Synonyms:	CT13; HAGE
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018665
ORF Size:	1944 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209792).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_018665.1</a>
RefSeq Size:	2318 bp
RefSeq ORF:	1947 bp
Locus ID:	55510
UniProt ID:	<a href="#">Q9NXZ2</a>
Cytogenetics:	6q13
MW:	72.7 kDa


[View online »](#)

**Gene Summary:**

The protein encoded by this gene is an ATP-dependent RNA helicase in the DEAD-box family and displays tumor-specific expression. [provided by RefSeq, Jul 2008]