

Product datasheet for **RC209807L4V**

ZNF706 (NM_016096) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ZNF706 (NM_016096) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ZNF706
Synonyms:	HSPC038; PNAS-106; PNAS-113
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_016096
ORF Size:	228 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209807).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016096.2
RefSeq Size:	2755 bp
RefSeq ORF:	231 bp
Locus ID:	51123
UniProt ID:	Q9Y5V0
Cytogenetics:	8q22.3
Domains:	zf-C2H2
MW:	8.5 kDa


[View online »](#)

Gene Summary:

Transcription repressor involved in the exit of embryonic stem cells (ESCs) from self-renewal.
Acts by repressing expression of KLF4.[UniProtKB/Swiss-Prot Function]