

Product datasheet for **RC209818L3V**

MSC (NM_005098) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MSC (NM_005098) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MSC
Synonyms:	ABF-1; ABF1; bHLHa22; MYOR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_005098
ORF Size:	618 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209818).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005098.3 , NP_005089.2
RefSeq Size:	2075 bp
RefSeq ORF:	621 bp
Locus ID:	9242
UniProt ID:	O60682
Cytogenetics:	8q13.3
Protein Families:	Druggable Genome, Transcription Factors
MW:	22.1 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene is a transcriptional repressor capable of binding an E-box element either as a homodimer or as a heterodimer with E2A in vitro. The encoded protein also forms heterodimers with E2A proteins in vivo. This protein is capable of inhibiting the transactivation capability of E47, an E2A protein, in mammalian cells. This gene is a downstream target of the B-cell receptor signal transduction pathway. [provided by RefSeq, Jul 2008]