

Product datasheet for **RC209851**

SRRD (NM_001013694) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SRRD (NM_001013694) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	SRRD
Synonyms:	HC/HCC; SRR1L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC209851 representing NM_001013694 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGCTGCGGCCGACGTGCGGCGCTGGAATCCTGGCAGGCGGCGGCTCCGCGAAGAGGCGCTCCGCGG
CTCGACGGCCGCGGCGGAGGGAGGCGGCCCGGGGAGAGAGGCGCGCCCGGGGAGAGAGGCGGC
GCCCCGGGGCCCCGAGGCGGAGTTCGAGTCTGACAGCGGAGTCGTGCTTCGTGCGCATCTGGGAGGCTGAG
AAGGACCTGTTTATCTCTGATTTCTGGAGTTCAGCACTAGAAACCATCAATAGATGTCTCACAAAACATC
TGAACAACCTGAAGGCCCTGTGGGACTCTTTCAGACATCTTTGAAACCTGCATCTTGACTCATTGCC
AGAGGAGTCAGATGTGGCACTGATTCTATCCCAAGAGAGATCTTGGTCACAGGAACCTGCCATTTGAAG
TGTGTGTGTTACGGCATTGGGAACCTTGCCACCTGCATCGTAGCTAGAAACCAAGTAACGTTTTTGTCTG
TTTTGTTGAAAAGTGCCAGATCCCAGAAGTCACTGTTGGGTATATGACCTCTGTTTAGCCAACCTTGA
AATTGAAGTCCTTAACACCCTTGGTGTGACTGTTCTCAGTGAGAACGAGGAAGGAAACGGAGTATTCGC
GGGAGCCTACCATCTTTACATGCTCCATTGTGGGACGGCCTTGTACAACATCTTTATGGAGTAAC
GGTCAGTAGATGCCCTTTCTAAGATGGTCATCATTGGGAACAGTTTCAAAGGACTTGAGGAGAGGTTGTT
GGCAAGGATTCTGCAGAAAAATTATCCCTACATTGCAAAGATTTTAAAGGACTGGAGGAGCTTGAGTTT
CCTCAGACTTCACAATACATGGACATATTTAATGATACCTCTGTCCACTGGTTCCCTGTGCAAAAGCTAG
AACAGCTCTCCATAGATATTTGGGAGTTTCGGGAAGAACCAGATTATCAGGACTGTGAGGACCTTGAAAT
CATCAGGAACAAGAGAGAAGATCCTTCTGCTACTGAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC209851 representing NM_001013694
 Red=Cloning site Green=Tags(s)


MAAAAAAALSWQAAAPRKRRSAARRPRRREAAPRGREAAPRGREAAPRGPEAFESDSGVVLRRIWEAE
 KDLFISDFWSSALETINRCLTKHLEQLKAPVGTLSDFGNLHLDLPEESDVATDSIPREILVTGTCHLK
 CVCYGIGNFATCIVARNQLTFLLLLLEKQIPRSHCWVDPLFSQLEIEVLNTLGVTVLSENEEGKRSIR
 GEPTIFYMLHCGTALYNNLLWSNWSVDALSKMVIIGNSFKGLEERLLARILQKNYPYIAKILKGLEEELEF
 PQTSQYMDIFNDTSVHWFPVQKLEQLSIDIWEFREEDYQDCEDLEIIRNKREDPSATD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8025_d03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001013694

ORF Size: 1017 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001013694.3](#)

RefSeq Size: 1310 bp

RefSeq ORF: 1020 bp

Locus ID: 402055


UniProt ID: [Q9UHQ3](#)

Cytogenetics: 22q12.1

MW: 38.4 kDa

Gene Summary: Plays a role in the regulation of heme biosynthesis and in the regulation of the expression of core clock genes.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC209851