

Product datasheet for **RC209863L2V**

PWWP domain containing protein 2B (PWWP2B) (NM_138499) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | PWWP domain containing protein 2B (PWWP2B) (NM_138499) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | PWWP domain containing protein 2B |
| Synonyms: | bA432J24.1; pp8607; PWWP2 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-mGFP (PS100071) |
| Tag: | mGFP |
| ACCN: | NM_138499 |
| ORF Size: | 1539 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC209863). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_138499.3 , NP_612508.2 |
| RefSeq Size: | 2626 bp |
| RefSeq ORF: | 1773 bp |
| Locus ID: | 170394 |
| UniProt ID: | Q6NUJ5 |
| Cytogenetics: | 10q26.3 |
| MW: | 55.7 kDa |


[View online »](#)

