

## Product datasheet for **RC209870L4V**

### DDX53 (NM\_182699) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	DDX53 (NM_182699) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DDX53
Synonyms:	CAGE; CT26
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_182699
ORF Size:	1893 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209870).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_182699.2</a>
RefSeq Size:	2151 bp
RefSeq ORF:	1896 bp
Locus ID:	168400
UniProt ID:	<a href="#">Q86TM3</a>
Cytogenetics:	Xp22.11
MW:	71.2 kDa


[View online »](#)

**Gene Summary:**

This intronless gene encodes a protein which contains several domains found in members of the DEAD-box helicase protein family. Other members of this protein family participate in ATP-dependent RNA unwinding. [provided by RefSeq, Sep 2011]