

## Product datasheet for **RC209878L2V**

### **SNX14 (NM\_153816) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	SNX14 (NM_153816) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SNX14
Synonyms:	RGS-PX2; SCAR20
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_153816
ORF Size:	2838 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209878).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_153816.2</a>
RefSeq Size:	3509 bp
RefSeq ORF:	2841 bp
Locus ID:	57231
UniProt ID:	<a href="#">Q9Y5W7</a>
Cytogenetics:	6q14.3
Protein Families:	Transmembrane
MW:	110 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the sorting nexin family. Members of this family have a phosphoinositide (PX) binding domain and are involved in intracellular trafficking. The encoded protein also contains a regulator of G protein signaling (RGS) domain. Regulator of G protein signaling family members are regulatory molecules that act as GTPase activating proteins for G alpha subunits of heterotrimeric G proteins. Alternate splicing results in transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2014]