

## Product datasheet for **RC209915L3V**

### ARHGAP17 (NM\_001006634) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ARHGAP17 (NM_001006634) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ARHGAP17
Synonyms:	MST066; MST110; MSTP038; MSTP066; MSTP110; NADRIN; PP367; PP4534; RICH-1; RICH1; WBP15
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001006634
ORF Size:	2643 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209915).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001006634.1</a>
RefSeq Size:	3526 bp
RefSeq ORF:	2646 bp
Locus ID:	55114
UniProt ID:	<a href="#">Q68EM7</a>
Cytogenetics:	16p12.1
Protein Families:	Druggable Genome



[View online »](#)

**MW:** 95.4 kDa

**Gene Summary:** RICH1 is a GTPase-activating protein (GAP). GAPs stimulate the intrinsic GTP hydrolysis of small G proteins, such as RHOA (MIM 165390), RAC1 (MIM 602048), and CDC42 (MIM 116952). [supplied by OMIM, Apr 2004]