

Product datasheet for **RC209917L4V**

ATF6 beta (ATF6B) (NM_004381) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ATF6 beta (ATF6B) (NM_004381) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ATF6 beta
Synonyms:	CREB-RP; CREBL1; G13
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_004381
ORF Size:	2109 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209917).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004381.4
RefSeq Size:	2622 bp
RefSeq ORF:	2112 bp
Locus ID:	1388
UniProt ID:	Q99941
Cytogenetics:	6p21.32
Protein Families:	Transcription Factors
MW:	76.7 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene is a transcription factor in the unfolded protein response (UPR) pathway during ER stress. Either as a homodimer or as a heterodimer with ATF6-alpha, the encoded protein binds to the ER stress response element, interacting with nuclear transcription factor Y to activate UPR target genes. The protein is normally found in the membrane of the endoplasmic reticulum; however, under ER stress, the N-terminal cytoplasmic domain is cleaved from the rest of the protein and translocates to the nucleus. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2008]