

Product datasheet for **RC209929A1V**

Human ILRUN (NM_024294) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ILRUN (NM_024294) AAV Particle
Tag:	Myc-DDK
Symbol:	ILRUN
Synonyms:	C6orf106; dj391O22.4; FP852
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC209929 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAGGGCATGGACGTAGACCTGGACCCGGAGCTGATGCAGAAGTTCAGCTGCCTGGGCACACCGACA
AGGACGTGCTCATCTCCGAGTTCAGAGGCTGCTCGGCTTCCAGCTCAATCCTGCCGTTGCGCCTTCTT
CCTGGACATGACCAACTGGAACCTACAAGCAGCAATTGGCGCCTATTATGACTTTGAGAGCCAAACATC
AGTGTGCCCTCTATGTCCTTTGTTGAAGATGTCACCATAGGAGAAGGGGAGTCAATACCTCCGGATACTC
AGTTTGTAAAAACATGGCGGATCCAGAATTCTGGGGCAGAGGCCTGGCCTCCAGGGGTTTGTCTTAAATA
TGTCGGGGGAGACCAATTTGGACATGTGAACATGGTGATGGTGAGATCGCTAGAGCCCCAAGAGATTGCA
GATGTCAGCGTCCAGATGTGCAGCCCCAGCAGAGCAGGAATGTATCAGGGACAGTGGCGGATGTGCACTG
CTACAGGACTCTACTATGGAGATGTATCTGGGTGATTCTCAGTGTGGAGGTGGGTGGACTTTTAGGAGT
AACGCAGCAGCTGTATCTTTTGAACGGAGTTCAACACACAGCCGCATCGTAAGGTAGAAGGAACTTC
AACCCTTTTGCCTCTCCCAAAAGAACCAGACAATCAGATGAAAACAACCTAAAAGACCCTGGGGGCTCCG
AGTTCGACTCGATCAGCAAAAACATGGGCTCCTGCTCCTGACACATGGGCTCCTGCTCCTGACCAAC
TGAGCAAGACCAGAATAGACTGTCACAGAACTCTGTAATCTGTCTCCAGCAGTCACGCAAAACAACCTTA
TCAGTAGTGACTTACAGTAAGGGGCTCCATGGGCTTACCCTTCGGCCAGTCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC209929 protein sequence Red=Cloning site Green=Tags(s) MEGMDVDLDPELMQKFSCLGTTDKDVLISEFQRLGFLNPAGCAFFLDMTNWNLQAAIGAYYDFESPNI SVPSMSFVEDVTIGEGESIPDQTQFVKTWRIQNSGAEAWPPGVCLKYVGGDQFGHVNMVMVRSLEPQEIA DYSVMCSPSRAGMYQGQWRMCTATGLYYGDVIWVILSVEVGGLLGVTQQLSSFETEFNTQPHRKVEGNF NPFASPQKNRQSDENNLKDPGGSEFDSISKNTWAPAPDTWAPAPDQTEQDQNRLSQNSVNLSPSSHANNL SVVTYSKGLHGPYPFGQS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_024294
ORF Size:	894 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_024294.2
RefSeq Size:	4433 bp
RefSeq ORF:	897 bp
Locus ID:	64771
UniProt ID:	Q9H6K1
Cytogenetics:	6p21.31
MW:	32.9 kDa