

Product datasheet for **RC209935L3V**

PXMP2 (NM_018663) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | PXMP2 (NM_018663) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | PXMP2 |
| Synonyms: | MPV17L3; PMP22 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_018663 |
| ORF Size: | 585 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC209935). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_018663.1 |
| RefSeq Size: | 977 bp |
| RefSeq ORF: | 588 bp |
| Locus ID: | 5827 |
| UniProt ID: | Q9NR77 |
| Cytogenetics: | 12q24.33 |
| Domains: | Mpv17_PMP22 |
| Protein Families: | Transmembrane |



[View online »](#)

MW: 22.3 kDa

Gene Summary: Seems to be involved in pore-forming activity and may contribute to the unspecific permeability of the peroxisomal membrane.[UniProtKB/Swiss-Prot Function]