

Product datasheet for RC209943A1V

Human ZFAND5 (NM_006007) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human ZFAND5 (NM_006007) AAV Particle
Tag: Myc-DDK
Symbol: ZFAND5
Synonyms: ZA20D2; ZFAND5A; ZNF216
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC209943 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCTCAGGAGACTAACCAGACCCCGGGGCCCATGCTGTGTAGCACAGGATGTGGCTTTTATGAAATC
 CTAGGACAAATGGAATGTGTTTCAGTTTGTCTACAAAGAATCTTCAGAGGCAGCAAAATAGTGCCAGAAT
 GAGCCCAATGGGACAGCTAGTGGTTCCAACAGTCCTACCTCAGATTCTGCATCTGTACAGAGAGCAGAC
 ACTAGCTTAACAACCTGTGAAGGTGCTGCTGGCAGCACATCTGAAAAATCAAGAAATGTGCCTGTGGCTG
 CCTTGCCTGTAACCTCAGCAAAATGACAGAAATGAGCATTTCAAGAGAGGACAAAATAACTACCCGAAAAC
 AGAGGTGTCAGAGCCAGTTGTCACTCAGCCAGTCCATCATGTTTCTCAGCCAGTACTTCTCAGAGTGAA
 GAAAAAGCTCCTGAATTGCCCAAAACAAAGAAAAACAGATGTTTCATGTGCAGAAAGAAAGTTGGTCTTA
 CAGGGTTTGACTGCCGATGTGGAAATTTGTTTGTGGACTTCACCGTTACTCTGACAAGCACAACGTGCC
 GTATGATTACAAAGCAGAAGCTGCAGCAAAATCAGAAAGAGAATCCAGTTGTTGTGGCTGAAAAAT
 CAGAGAATA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC209943 protein sequence
 Red=Cloning site Green=Tags(s)

MAQETNQTPGMLCSTGCGFYGNPRTNGMCSVCYKEHLQRQQNSGRMSPMGTASGSNSPTSDSASVQRAD
 TSLNCEGAAGSTSEKSRNPVAALPVTQQMTESISREDKITTPKTEVSEPVVTQSPSVSQPSTSQSE
 EKAPELPKPKKNRCFMCRKKVGLTGFDRCGNLFCGLHRYSDKHNCPTYDYKAEAAKIRKENPVVVAEKI
 QRI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_006007
ORF Size:	639 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_006007.1
RefSeq Size:	5747 bp
RefSeq ORF:	642 bp
Locus ID:	7763
UniProt ID:	O76080
Cytogenetics:	9q21.13
MW:	23.1 kDa