

Product datasheet for RC209945A1V

Human P311 (NREP) (NM_004772) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human P311 (NREP) (NM_004772) AAV Particle
Tag:	Myc-DDK
Symbol:	P311
Synonyms:	C5orf13; D4S114; P311; PRO1873; PTZ17; SEZ17
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	<p>>RC209945 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)</p> <p>TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC</p> <p>ATGGTTTATTACCCAGAACTCTTTGTCTGGTCAAGAACCAATTTCCAAACAAGGACATGGAGGGAA GGCTTCCTAAGGGAAGACTTCCTGTCCCAAAGGAAGTGAACCGCAAGAAGACGATGAGACAAACGCTGC CTCCCTGACTCCACTGGGCAGCAGTGAAGTCCGCTCCCAAGAATCAGTTACCTCCACTTTTT</p> <p>ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT ACAAGGATGACGACGATAAGGTTTAA</p>
Protein Sequence:	<p>>RC209945 protein sequence Red=Cloning site Green=Tags(s)</p> <p>MVYYPELFWVVSQEPFPNKDMEGRLLPKGRLLPVPKEVNRKKNDETNAASLTPLGSELRSRPRISYLHFF</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Species:	Human
Serotype:	AAV-2
ACCN:	NM_004772
ORF Size:	204 bp
Buffer:	PBS with 0.001% Pluronic F68



[View online »](#)

Stability: AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

RefSeq: [NM_004772.1](#)

RefSeq Size: 2026 bp

RefSeq ORF: 207 bp

Locus ID: 9315

UniProt ID: [Q16612](#)

Cytogenetics: 5q22.1

MW: 7.9 kDa