

## Product datasheet for **RC209950L3V**

### Prostaglandin I2 Receptor (PTGIR) (NM\_000960) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Prostaglandin I2 Receptor (PTGIR) (NM_000960) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Prostaglandin I2 Receptor
Synonyms:	IP; PRIPR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_000960
ORF Size:	1158 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209950).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_000960.3</a>
RefSeq Size:	2103 bp
RefSeq ORF:	1161 bp
Locus ID:	5739
UniProt ID:	<a href="#">P43119</a>
Cytogenetics:	19q13.32
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Neuroactive ligand-receptor interaction, Vascular smooth muscle contraction



[View online »](#)

**MW:** 41 kDa

**Gene Summary:** The protein encoded by this gene is a member of the G-protein coupled receptor family 1 and has been shown to be a receptor for prostacyclin. Prostacyclin, the major product of cyclooxygenase in macrovascular endothelium, elicits a potent vasodilation and inhibition of platelet aggregation through binding to this receptor. [provided by RefSeq, Jul 2008]