

Product datasheet for **RC209971L1V**

ZADH2 (NM_175907) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | ZADH2 (NM_175907) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | ZADH2 |
| Synonyms: | PRG-3 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_175907 |
| ORF Size: | 1131 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC209971). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_175907.3 , NP_787103.1 |
| RefSeq Size: | 5495 bp |
| RefSeq ORF: | 1134 bp |
| Locus ID: | 284273 |
| UniProt ID: | Q8N4Q0 |
| Cytogenetics: | 18q22.3 |
| Protein Families: | Druggable Genome |
| MW: | 40.1 kDa |



[View online »](#)

Gene Summary:

Functions as 15-oxo-prostaglandin 13-reductase and acts on 15-keto-PGE1, 15-keto-PGE2, 15-keto-PGE1-alpha and 15-keto-PGE2-alpha with highest efficiency towards 15-keto-PGE2-alpha. Overexpression represses transcriptional activity of PPARG and inhibits adipocyte differentiation.[UniProtKB/Swiss-Prot Function]