

Product datasheet for **RC209973A1V**

Human IL17F (NM_052872) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human IL17F (NM_052872) AAV Particle
Tag: Myc-DDK
Symbol: IL17F
Synonyms: CANDF6; IL-17F; ML-1; ML1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC209973 representing NM_052872
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGACAGTGAAGACCCTGCATGGCCAGCCATGGTCAAGTACTTGCTGCTGTCGATATTGGGGCTTGCTT
TTCTGAGTGAGGCGGCAGCTCGGAAAATCCCAAAGTAGGACATACTTTTTCCAAAAGCCTGAGAGTTG
CCCGCCTGTGCCAGGAGGTAGTATGAAGCTTGACATTGGCATCATCAATGAAAACCGCGCTTTCCATG
TCACGTAACATCGAGAGCCGCTCCACCTCCCCCTGGAATTACACTGTCACTTGGGACCCCAACCGGTACC
CCTCGGAAGTTGTACAGGCCAGTGTAGGAACCTGGGCTGCATCAATGCTCAAGGAAAGGAAGACATCTC
CATGAATTCCGTTCCCATCCAGCAAGAGACCCTGGTCCGCGAGGAAGCACCAAGGCTGCTCTGTTTCT
TTCCAGTTGGAGAAGGTGCTGGTGAAGTGTGGCTGCACCTGCGTCACCCCTGTCATCCACCGTGTGCAG

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC209973 representing NM_052872
Red=Cloning site **Green**=Tags(s)

MTVKTLHGPMVKYLLLSILGLAFLSEAAARKIPKVGHTFFQKPESCPPVPGGSMKLDIGIINENQRVSM
SRNIESRSTSPWNYTVTWDPNRYSPSEVVQAQCRNLGCINAQKEDISMNSVPIQQETLVVRKHQGSVS
FQLEKVLVTGCTCVTPVIHRVQ

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online »](#)

ACCN:	NM_052872
ORF Size:	489 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_052872.3</u>
RefSeq Size:	808 bp
RefSeq ORF:	492 bp
Locus ID:	112744
UniProt ID:	<u>Q96PD4</u>
Cytogenetics:	6p12.2
MW:	14.9 kDa