

Product datasheet for **RC209991L3V**

HOXB13 (NM_006361) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	HOXB13 (NM_006361) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HOXB13
Synonyms:	HPC9; PSGD
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_006361
ORF Size:	852 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC209991).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006361.5
RefSeq Size:	3047 bp
RefSeq ORF:	855 bp
Locus ID:	10481
UniProt ID:	Q92826
Cytogenetics:	17q21.32
Domains:	homeobox
Protein Families:	Transcription Factors



[View online »](#)

MW: 30.7 kDa

Gene Summary: This gene encodes a transcription factor that belongs to the homeobox gene family. Genes of this family are highly conserved among vertebrates and essential for vertebrate embryonic development. This gene has been implicated to play a role in fetal skin development and cutaneous regeneration. In mice, a similar gene was shown to exhibit temporal and spatial colinearity in the main body axis of the embryo, but was not expressed in the secondary axes, which suggests functions in body patterning along the axis. This gene and other HOXB genes form a gene cluster at chromosome the 17q21-22 region. [provided by RefSeq, Jul 2008]