

Product datasheet for **RC210003A1V**

Human CCR4 (NM_005508) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human CCR4 (NM_005508) AAV Particle
Tag:	Myc-DDK
Symbol:	CCR4
Synonyms:	CC-CKR-4; CD194; ChemR13; CKR4; CMKBR4; HGCN:14099; K5-5
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC210003 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAACCCACGGATATAGCAGACACCACCTCGATGAAAGCATATACAGCAATTACTATCTGTATGAAAGTATCCCAAGCCTTGCACCAAGAAGGCATCAAGGCATTTGGGGAGCTCTTCTGCCCCCACTGTATTCCTTGGTTTTGTATTTGGTCTGCTTGGAAATCTGTGGTGGTTCTGGTCTGTCAAATACAAGCGGCTCAGGTCCATGACTGATGTGTACCTGCTCAACCTTGCCATCTCGGATCTGCTCTTCGTGTTTTCCCTCCCTTTTGGGGCTACTATGCAGCAGACCAAGTGGGTTTTTGGGCTAGGTCTGTGCAAGATGATTTCTGGATGTACTTGGTGGGCTTTTACAGTGGCATATTCTTTGTCATGCTCATGAGCATTGATAGATACCTGGCAATTGTGCACGCGGTGTTTTCTTGAGGGCAAGGACCTTGACTTATGGGGTCATCACCAGTTTGGCTACATGGTCAGTGGCTGTGTTCCCTCCCTTCTGGCTTTCTGTTGAGCACTTGTATGAGCGCAACCATACCTACTGCAAAACCAAGTACTCTCTCAACTCCACGACGTGGAAGTTCTCAGCTCCCTGGAAATCAACATTCTCGGATTGGTGATCCCCTTAGGGATCATGCTGTTTGTACTCCATGATCATCAGGACCTTGACGATTGAAAAATGAGAAGAAGAACAAGCGGTGAAGATGATCTTGCCGTGGTGGTCTCTTCTTGGGTTCTGGACACCTTACAACATAGTGCTCTTCTAGAGACCCTGGTGGAGCTAGAAGTCTCAGGACTGCACCTTTGAAAGTACTTGGACTATGCCATCCAGGCCACAGAACTCTGGCTTTTGTCACTGCTGCCTTAATCCCATCATCTACTTTTTCTGGGGGAGAAATTCGCAAGTACATCCTACAGCTCTTCAAACCTGCAGGGGCCTTTTGTGCTCTGCCAATACTGTGGGCTCCTCAAATTTACTCTGCTGACACCCCACTCATCTTACACGCAGTCCACCATGGATCATGATCTCCATGATGCTCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC210003 protein sequence Red=Cloning site Green=Tags(s) MNPTDIADTTLDESIYSNYLYESIPKPCTKEGIKAFGELFLPPLYSLVFVFGLLGNSVVVLVLFKYKRL RSMTDVYLLNLAISDLLFVFSLPFWGYAADQWVFGGLCKMISWMLVGFYSGIFFVMLMSIDRYLAIV HAVFSLRARTLYGVITSLATWSVAVFASLPGFLFSTCYTERNHTYCKTKYSLNSTTWKVLSSLEINILG LVIPLGIMLFCYSMIIRTLQHCKNEKKNKAVKMIFAVVVLFLGFWTPYNIIVLFLETLELEVLQDCTFER YLDYAIQATETLAFVHCCLNPPIIYFFLGEKFRKYILQLFKTCRGLFVLCQYCGLLQIYSADTPSSSYTQS TMDHDLHDAL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_005508
ORF Size:	1080 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_005508.4
RefSeq Size:	1657 bp
RefSeq ORF:	1083 bp
Locus ID:	1233
UniProt ID:	P51679
Cytogenetics:	3p22.3
MW:	41.4 kDa