

Product datasheet for RC210004L1V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

CD40L (CD40LG) (NM_000074) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: CD40L (CD40LG) (NM_000074) Human Tagged ORF Clone Lentiviral Particle

Symbol: CD40L

Synonyms: CD40L; CD154; gp39; hCD40L; HIGM1; IGM; IMD3; T-BAM; TNFSF5; TRAP

Mammalian Cell

Selection:

None

Vector: pLenti-C-Myc-DDK (PS100064)

Tag: Myc-DDK
ACCN: NM 000074

ORF Size: 783 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC210004).

Sequence:

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This

naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeq: <u>NM 000074.2</u>

RefSeq Size: 1834 bp
RefSeq ORF: 786 bp
Locus ID: 959
UniProt ID: P29965
Cytogenetics: Xq26.3

Domains: TNF

Protein Families: Druggable Genome, Secreted Protein, Transmembrane





CD40L (CD40LG) (NM_000074) Human Tagged ORF Clone Lentiviral Particle - RC210004L1V

Protein Pathways: Allograft rejection, Asthma, Autoimmune thyroid disease, Cell adhesion molecules (CAMs),

Cytokine-cytokine receptor interaction, Primary immunodeficiency, Systemic lupus

erythematosus, T cell receptor signaling pathway, Viral myocarditis

MW: 29.3 kDa

Gene Summary: The protein encoded by this gene is expressed on the surface of T cells. It regulates B cell

function by engaging CD40 on the B cell surface. A defect in this gene results in an inability to undergo immunoglobulin class switch and is associated with hyper-lgM syndrome. [provided

by RefSeq, Jul 2008]