

Product datasheet for **RC210032A1V**

Human Tetraspanin 9 (TSPAN9) (NM_006675) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Tetraspanin 9 (TSPAN9) (NM_006675) AAV Particle
Tag:	Myc-DDK
Symbol:	Tetraspanin 9
Synonyms:	NET-5; NET5; PP1057
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC210032 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCAGGGGCTGCCTCTGCTGCTTGAAGTACATGATGTTCTCTTCAATTTGATATTCTGGCTCTGTG
GCTGTGGGCTGCTGGGAGTGGGCATCTGGCTCTCCGTGTCCCAAGGCACTTTGCCACCTTCTCCCCAG
CTTCCCTTCGTTGTCTGCAGCCAACCTGGTCATTGCCATAGGCACCATTTGTCATGGTGACGGGCTTCTC
GGCTGCCTGGGGGCCATCAAGGAAAACAAGTGCCTCCTCCTCAGCTTTTTCATCGTCCTGTTGGTCATCC
TCCTAGCAGAGCTGATCTTACTCATCCTCTTCTTTGTCTACATGGACAAGGTGAACGAGAACGCCAAGAA
GGACCTGAAGGAAGGCTGCTGTGTACACACCCGAGAACAACGTGGGGCTGAAGAACGCCTGGAACATC
ATCCAGGCTGAGATGCGATGCTGTGGTGTCACTGACTACACAGACTGGTACCCAGTGCTGGGGGAGAACA
CGGTTCCCGACCGCTGCTGCATGGAGAACTCCAGGGCTGCGGGCGCAACGCCACCACGCCTTTGTGGAG
AACGGGCTGCTATGAAAAGGTGAAGATGTGGTTCGATGACAATAAGCACGTGCTGGGCACGGTGGGGATG
TGCATCCTCATCATGCAGATCCTGGGCATGGCCTTCTCCATGACCCTCTCCAGCACATCCACCGGACTG
GTAAGAAGTACGACGCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC210032 protein sequence Red=Cloning site Green=Tags(s) MARGCLCCLKYMMFLFNLIFWLCGCLLGVGIWLSVSQGNFATFSPSFPSLSAANLVIAIGTIVMTGFL GCLGAIKENKCLLLSFFIVLLVILLAELILLILFFVYMDKVNENAKDLKEGLLLYHTENNVGLKNAWNI IQAEMRCCGVTDYTDWYPVLGENTVPDRCCMENSQGCGRNATTPLWRTGCEYKVKMWFFDDNKHVLGTVGM CILIMQILGMAFSMTLFQHIHRTGKKYDA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_006675
ORF Size:	717 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_006675.3</u>
RefSeq Size:	4351 bp
RefSeq ORF:	720 bp
Locus ID:	10867
UniProt ID:	<u>O75954</u>
Cytogenetics:	12p13.33-p13.32
MW:	26.8 kDa