

## Product datasheet for **RC210036L3V**

### SOX7 (NM\_031439) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SOX7 (NM_031439) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SOX7
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031439
ORF Size:	1164 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210036).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031439.2</a>
RefSeq Size:	3278 bp
RefSeq ORF:	1167 bp
Locus ID:	83595
UniProt ID:	<a href="#">Q9BT81</a>
Cytogenetics:	8p23.1
Domains:	HMG
Protein Families:	ES Cell Differentiation/IPS, Induced pluripotent stem cells, Transcription Factors
MW:	42.2 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional regulator after forming a protein complex with other proteins. The protein may play a role in tumorigenesis. A similar protein in mice is involved in the regulation of the wingless-type MMTV integration site family (Wnt) pathway. [provided by RefSeq, Jul 2008]