

Product datasheet for **RC210050L4V**

Laminin 5 (LAMB3) (NM_001017402) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Laminin 5 (LAMB3) (NM_001017402) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | LAMB3 |
| Synonyms: | A11A; BM600-125KDA; LAM5; LAMNB1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001017402 |
| ORF Size: | 3516 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC210050). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001017402.1 , NP_001017402.1 |
| RefSeq Size: | 4339 bp |
| RefSeq ORF: | 3519 bp |
| Locus ID: | 3914 |
| UniProt ID: | Q13751 |
| Cytogenetics: | 1q32.2 |
| Protein Families: | Druggable Genome, Secreted Protein |
| Protein Pathways: | ECM-receptor interaction, Focal adhesion, Pathways in cancer, Small cell lung cancer |



[View online »](#)

MW: 129.6 kDa

Gene Summary: The product encoded by this gene is a laminin that belongs to a family of basement membrane proteins. This protein is a beta subunit laminin, which together with an alpha and a gamma subunit, forms laminin-5. Mutations in this gene cause epidermolysis bullosa junctional Herlitz type, and generalized atrophic benign epidermolysis bullosa, diseases that are characterized by blistering of the skin. Multiple alternatively spliced transcript variants that encode the same protein have been found for this gene. [provided by RefSeq, Jul 2008]