

Product datasheet for **RC210080L3V**

RISC (SCPEP1) (NM_021626) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	RISC (SCPEP1) (NM_021626) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RISC
Synonyms:	HSCP1; RISC
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_021626
ORF Size:	1356 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC210080).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_021626.1
RefSeq Size:	1952 bp
RefSeq ORF:	1359 bp
Locus ID:	59342
UniProt ID:	Q9HB40
Cytogenetics:	17q22
Domains:	serine_carbpept
Protein Families:	Druggable Genome, Protease, Secreted Protein



[View online »](#)

MW: 50.8 kDa

Gene Summary: May be involved in vascular wall and kidney homeostasis.[UniProtKB/Swiss-Prot Function]