

## Product datasheet for **RC210100**

### TARP (NM\_001003799) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** TARP (NM\_001003799) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** TARP  
**Synonyms:** CD3G; TCRG; TCRGC1; TCRGC2; TCRGV  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >RC210100 representing NM\_001003799  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGCAGATGTTTCCCCAAGCCCACTATTTTTCTTCTCAATTGCTGAAACAAAGCTCCAGAAGGCTGG  
 AACATACCTTTGTCTTCTTGAGAAATTTTCCCTGATGTTATTAAGATACATTGGCAAGAAAAGAAGAGC  
 AACACGATTCTGGGATCCAGGAGGGGAACACCA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC210100 representing NM\_001003799  
 Red=Cloning site Green=Tags(s)

MQMFPPSPLFFFLQLLKQSSRRLEHTFVFLRNFSMLLLRYIGKRRRATRFWDPRRGTP


**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001003799  
 ORF Size: 174 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001003799.1</a></u> , <u><a href="#">NP_001003799.1</a></u>
<b>RefSeq Size:</b>	1027 bp
<b>RefSeq ORF:</b>	177 bp
<b>Locus ID:</b>	445347
<b>Cytogenetics:</b>	7p14.1
<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	7.2 kDa
<b>Gene Summary:</b>	In some non-lymphoid tissues, the unrearranged T cell receptor gamma (TRG@) locus is expressed. The resulting transcript contains a subset of the TRG@ gene segments and is shorter than TRG@ transcripts expressed in lymphoid tissues. This RefSeq record represents the unrearranged TRG@ locus transcript; the complete TRG@ locus is represented by the genomic RefSeq NG_001336. The transcript represented by this RefSeq has two open reading frames (ORFs) that encode different proteins. The downstream ORF is in the same frame as TRG@ and its protein product is similar to TRG@ proteins. The upstream ORF uses a different reading frame and encodes a novel protein. [provided by RefSeq, Jul 2008]