

Product datasheet for **RC210118A1V**

Human Caveolin 3 (CAV3) (NM_001234) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human Caveolin 3 (CAV3) (NM_001234) AAV Particle
Tag: Myc-DDK
Symbol: Caveolin 3
Synonyms: LGMD1C; LQT9; MPDT; RMD2; VIP-21; VIP21
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC210118 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGATGGCAGAAGAGCACACAGATCTCGAGGCCAGATCGTCAAGGATATCCACTGCAAGGAGATTGACC
 TGGTGAACCGAGACCCCAAGAACATTAAATGAGGACATAGTCAAGGTGGATTTGAAGACGTGATCGCAGA
 GCCTGTGGGCACCTACAGCTTTGACGGCGTGTGGAAGGTGAGCTACACCACCTTCACTGTCTCCAAGTAC
 TGGTGTACCGTCTGTTGTCCACGCTGCTGGGCGTCCCACTGGCCCTGCTCTGGGGCTTCTGTTCCGCT
 GCATCTCCTTCTGCCACATCTGGGCGGTGGTGCCATGCATTAAGAGCTACCTGATCGAGATCCAGTGCAT
 CAGCCACATCTACTCACTCTGCATCCGACCTTCTGCAACCCACTCTTCGCGGCCCTGGGCCAGGTCTGC
 AGCAGCATCAAGGTGGTGTGCTGCGGAAGGAGGTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC210118 protein sequence
 Red=Cloning site Green=Tags(s)

MMAEEHTDLEAQIVKDIHCKEIDLVRDPKINEDIVKVDFEDVIAEPVGTYSFDGVWVSYTFTVSKY
 WCYRLLSTLLGVPLALLWGFLFACISFCHIWAIVPCIKSYLIEIQCISHIYSLCIRTFNPLFAALGQVC
 SSIKVVLRKEV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2


[View online »](#)

ACCN:	NM_001234
ORF Size:	453 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001234.3</u>
RefSeq Size:	1339 bp
RefSeq ORF:	456 bp
Locus ID:	859
UniProt ID:	<u>P56539</u>
Cytogenetics:	3p25.3
MW:	17.3 kDa