

Product datasheet for **RC210136A1V**

Human Kallikrein 13 (KLK13) (NM_015596) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Kallikrein 13 (KLK13) (NM_015596) AAV Particle
Tag:	Myc-DDK
Symbol:	Kallikrein 13
Synonyms:	KLK-L4; KLKL4
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC210136 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTGGCCCTGGCCCTAGTGATCGCCTCCCTGACCTTGGCCTTGTCAGGAGGTGTCTCCAGGAGTCTT
CCAAGTTCTCAACACCAATGGGACCAAGTGGTTTCTCCAGGTGGCTACACCTGCTTCCCCACTCTCA
GCCCTGGCAGGCTGCCCTACTAGTGCAAGGGCGGCTACTCTGTGGGGAGTCTGGTCCACCCCAATGG
GTCCTCACTGCCGCACACTGTCTAAAGGAGGGGCTCAAAGTTTACCTAGGCAAGCACGCCCTAGGGCGTG
TGAAGCTGGTGAGCAGGTGAGGGAAGTTGTCCACTCTATCCCCACCCTGAATACCGGAGAAGCCCCAC
CCACCTGAACCACGACCATGACATCATGCTTCTGGAGCTGCAGTCCCCGGTCCAGCTCACAGGTACATC
CAAACCTGCCCCCTTTCCACAACAACCGCCTAACCCCTGGCACCACCTGTCGGGTGTCTGGCTGGGGCA
CCACCACGAGCCCCAGGTGAATTACCCAAAACCTCTACAATGTGCCAACATCCAACCTCGCTCAGATGA
GGAGTGTGCTCAAGTCTACCCAGGAAAGATCACTGACAACATGTTGTGTGCCGGCACAAAGAGGGTGGC
AAAGACTCCTGTGAGGGTGACTCTGGGGCCCCCTGGTCTGTAAACAGAACTGTATGGCATCGTCTCT
GGGAGACTTCCCATGTGGGCAACCTGACCGGCCTGGTGTCTACACCCGTGTCTAAGATACGCTCTGTG
GATCCGTGAAACAATCCGAAAAATGAAACCCAGCAGCAAAAATGTTGAAGGGCCACAA

ACGCGTACGCGGCGGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC210136 protein sequence Red=Cloning site Green=Tags(s) MWPLALVIASLTALSGGVSQESSKVLNTNGTSGFLPGGYTCFPHSQPWQAALLVQGRLLCGGVLVHPKW VLTAACHLKEGLKVYLKGHALGRVEAGEQVREVVHSIPHPEYRRSPTHLNHDHIMLLELQSPVQLTGYI QTLPLSHNNRLTPGTTTCRVSGWGTSTSPQVNYPKTLQCANIQLRSDEECRQVYPGKITDNMLCAGTKEGG KDSCEGDSGGPLVCNRTLYGIVSWGDFPCGQPDPRPGVYTRVSRVVLWIRETIRKYETQQQKWLKGPQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_015596
ORF Size:	831 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_015596.1</u>
RefSeq Size:	1258 bp
RefSeq ORF:	834 bp
Locus ID:	26085
UniProt ID:	<u>Q9UKR3</u>
Cytogenetics:	19q13.41
MW:	30.6 kDa